Liners & Sleeves

Regulatory notes for professionals

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ENGLISH

*The word « product » used in this document applies to both Your® Liners and Your® Sleeves.

1. Troubleshooting

MotionTech's ordering platform allows you to report issues and file warranty claims. It also contains instructions and guidance on potential issues you might face. <https://motiontech.ch/yourlinerorderingplatform>

Liner issues

Problem	Likely cause	Possible solution
Seam is wearing or unraveling	Rough socket edge and/or rubbing against the socket.	Apply silicone glue to the seam
Fabric comes off silicone at proximal edge	Liner is trimmed too short and rubs against the socket edge	Liner needs to be changed and trimmed longer
Fabric shows signs of wear along socket trim line	The liner rubs against the trim line of the socket	Smooth out the socket edge and cover it with protective layer.
Fabric shows signs of wear	Socket has become loose in that area	Tighten the socket to avoid friction between the liner and the socket
Cracks in the silicone	Liner exposed to excessive pressure in this area	Relieve the pressure in this area by reworking the socket fit If the problem persists, using a firmer silicone may help to solve the problem.
Excessive sweating	Skin is not yet used to the closed environment that a liner imposes	Excessive sweating should reduce after a few weeks of using Your® Liner.
Excessive sweating appears although not experienced at first	Some air is getting into the liner	This may be due to limb shrinkage. Check the socket and liner fit. If the liner fit is loose, order a new Your® Liner with updated scan. If the socket fit is loose, tighten up the socket.
	Liner presents a hole	Your® Liner needs to be renewed
Top edge of the liner is rolling down the leg	Distal soft tissues not well accounted for in 3D model	Compressible soft tissues might offset how the liner is on the limb compared to how it was intended. Contact MotionTech service team.
	Unadapted proximal compressions	The 3D model we made is too tight at proximal edge of liner. Contact MotionTech service team.
	Limb volume increased	The limb may have increased in volume. A new Your® Liner shall be ordered with updated scan.

Problem	Likely cause	Possible solution
Patient presents irritation all over the residual limb	Residual limb has shrunk and the liner or socket fit became loose	The socket fit shall be tightened. If this is not efficient, a new Your® Liner shall be ordered with updated scan for the new volume and shape of your patient's limb. To ensure that the irritation does not come from allergic reaction, apply a piece of silicone from a Your® Liner in another area of the body and check for skin reaction.
	Use of an irritating soap lotion, cream, or other skin conditioner	Ask your patient if he started using a new cleaning product, lotion, cream, powder or other skin conditioner. Hydrocarbon oils, animal fats and oils can often produce allergic reactions. Antibacterial soap has also been linked to skin irritations with some amputees.
	Hygiene routine not properly implemented	Review the hygiene instructions with your patient, including cleaning, rinsing and disinfection procedures.
	Patient does not follow the proper rolling procedure and pulls the liner up on the the skin instead.	Check your patient's donning technique. Reduced irritations have also been obtained by cutting the proximal edge of the liner in a wave pattern.
Open wounds or other skin injuries	Many possible causes. Socket fit and other specific patients conditions are most likely.	The specific cause need to be singled out. During healing, Your® Liner should not be worn if the condition is too severe. If the injury is minor, bandages should be applied to cover the wound and prevent direct contact with Your® Liner.
Residual limb and liner rotate inside the socket	Residual limb is highly cylindrical and socket is too loose	Tighten socket fit. Re-order a Your® Liner and ask to add anti-rotation silicone rings.

2. Contraindications

Contraindications for all Your® Products

Although silicone Liners and Sleeves have been used successfully with thousands of patients, MotionTech does not recommend using Your® Liners/Sleeves with the following types of patients:

- any patient who is either not capable or not willing to comply with the cleaning and hygiene procedures described in these Instructions for Use
- any patient who lacks the hand function to roll the product successfully up onto the residual limb, and who does not have a caregiver to assist with this process
- any patient who lacks the ability to orient the product properly on the limb, for example due to limited eyesight
- any patient who has a skin condition which could be aggravated by covering the limb with a non-porous product, or who is known for developing allergic reactions to silicone

Contraindications for Your® Suction Liners

Your® Suction Liners are intended to work safely when used with a dedicated prosthetic socket equipped with a standard expulsion valve. A Your® Suction Liner must also be used with a supplementary suspension system independent from the liner, in case the patient presents a short residual limb, substantial daily volume fluctuations, or if sufficient suspension can not be achieved without such supplementary suspension system. A

Your® Suction Liner can not be provided without supplementary suspension method to a patient who presents a short residual limb or substantial residual limb volume fluctuations during the day, due to the risks of loss of suspension.

Contraindications for Your® Femoral Pin-lock Liners

Your® Femoral Pin-lock Liners shall not be provided to patients with a very short residual limb (approximately 7cm (3") or shorter as measured from the perineum to the distal end along the medial side), unless a supplementary suspension system independent from the liner is used (for example a Silesian belt or a pelvic band), due to the risks of loss of suspension.

Contraindications for Your® Sleeve from Measurements

A Your® Sleeve from Measurements is designed to function safely when positioned to provide secure suspension of the prosthesis to the limb. Because of the risk of air leakage which could induce a loss of suspension, a Your® Sleeve from Measurements should not be supplied to a patient who presents any of the under contraindications:

- . Invaginated scars or longitudinal concavities able to form an air channel which would lead to improper hold of the vacuum
- . A too short residual limb (i.e. less than 10cm)
- . A limb with a non-rectilinear central axis (e.g. strong varus/valgus, certain types of congenital limbs).

3. Suspension type & compatible components

Suspension type & components compatible with Your® Tibial Cushion Liners

A Your® Tibial Cushion Liner is designed to be used with a suspension sleeve. Even if you decide to add some anatomical suspension in your socket design - which is not recommended, as Your® Tibial Cushion Liners work best with total weight bearing sockets -, you MUST provide a suspension sleeve to your patient in order for the suspension to be efficient. Not using a suspension sleeve with a Your® Tibial Cushion Liners may result in slippage of the prosthesis on the liner.

It is your responsibility, as ordering prosthetist, to ensure that the suspension function provided by the chosen suspension sleeve is adequate and sufficient.

A Your® Tibial Cushion Liner does not hold any weight bearing or suspension function, and its strength is therefore not safety-related. Your® Tibial Cushion Liners have therefore not been subjected to strength tests. For further indications on maximum prosthesis weight and/or user weight, refer to the indications provided by the manufacturer of the components involved in the weight bearing and suspension functions.

Suspension type & components compatible with Your® Suction Liners

A Your® Suction Liner is designed to be used with a dedicated socket equipped with a standard expulsion valve. Before choosing this type of liner for a patient, observe the contraindications presented in the present manual. It is mandatory that you ensure that the suspension strength obtained through the vacuum suspension is sufficient prior to providing the Your® Suction Liner to your patient. If sufficient suspension can not be reached through socket adaptations, you must provide him with a supplementary suspension system independent from the liner.

Extensive efforts have been implemented to ensure that the suspension strength of Your® Suction Liners is as high as possible. However, as Your® Suction Liners are tailor-made, classical strength tests are not sufficient to ensure a flawless mechanical strength for each designed product. Furthermore, the strength of the suspension depends to a great extent on the fit with the prosthetic socket. This is why the suspension function of Your® Suction Liners must be validated by the certified prosthetist prior to provision of such a device to the patient.

Your® Suction Liners have therefore only been subjected to shear fatigue tests to assess the strength of the silicone rings. If any sign of excessive movement or product wear is detected on or around the silicone rings of the liner at any time, your patient should be instructed to immediately discontinue use of the device and consult you.

For further indications on maximum prosthesis weight and/or user weight, refer to the indications provided by the manufacturer of the components involved in the weight bearing and suspension functions.

Suspension type & components compatible with Your® Pin-lock Liners

A Your® Pin-lock Liner is designed to be used with a distal attachment pin. MotionTech provides two options regarding compatible pin-lock systems: a M6 and a M10 thread. It is your responsibility, as ordering prosthetist, to ensure that the selected pin-lock system to be assembled on the prosthesis is compatible with the chosen connector (M6 or M10).

Apply Loctite 410/411 or equivalent compound to the threads of the pin and torque to 4 Nm (3 ft.lbs). It is not recommended to use

a pin without a flange with Your® Pin-lock Liners. If you do, follow the same process but torque only to 1.5 Nm (11 ft.lbs).

Once Loctite the pin has been fixated with Loctite, it must not be removed from the liner. Trying to remove the pin from the liner might result in weakening the umbrella, which might lead to failure of the umbrella during use, exposing your patient to the risk of falling.

Extensive efforts have been implemented to ensure that the suspension strength of Your® Pin-lock liners is as high as possible. Your® Pin-lock Liners have been tested under axial tensile load to assess their proof and ultimate strength under such loading conditions. However, as the shape of the locking system of Your® Pin-lock Liners are tailored to each patient, classical strength tests are not sufficient to ensure a flawless mechanical strength for each designed product. These conducted tests provided indications that a Your® Pin-lock Liner may be used with a prosthesis which weighs up to 10kg.

If any sign of excessive movement or product wear is detected on or around the locking part of the liner at any time, your patient should be instructed to immediately discontinue use of the device and consult you.

For further indications on maximum prosthesis weight and/or user weight, refer to the indications provided by the manufacturer of the components involved in the weight bearing and suspension functions.

Components compatible with Your® Sleeve

A Your® Sleeve is designed to be used with a dedicated socket equipped with a standard expulsion valve.

There is no indications on maximum prosthesis weight, when using a Your® Sleeve as the secure suspension function depends on the ability of the sleeve to act as an airtight seal and not on its ability to carry weight.

4. Product reception & fitting

Product size check

Despite all the quality checks we implement throughout our production, MotionTech products remain custom-made devices. It is therefore mandatory that you be present the first time the product is donned onto your patient's limb to ensure its suitability before providing it to your patient.

Pay attention to the observed tension of the product once donned and make sure your patient does not feel unusual compression on his/her limb. Leave the product on your patient's limb for at least 5 minutes to ensure that the compression is neither too high nor too low. If there is a doubt on the compressive behaviour of the product, contact MotionTech's customer service prior to letting your patient go with the product. If air pockets or wrinkles persist to appear after five attempts of carefully donning and aligning the product, then the product should be deemed not adapted, and you should contact MotionTech' customer service for further indications. Similarly, if the product is too tight, too loose or too short, you should contact MotionTech' customer service to assess the situation.

Checking the suspension of a Your® Sleeve

When using a Your® Sleeve, please also ensure that the suspension it brings to the prosthesis is safe and adequate. First test by simply pulling on the prosthesis while the patient is sitting. Once you are sure the suspension is safe, let your patient test the suspension by walking gradually longer walks in a controlled environment, ensuring you avoid the risks of falling.

If in doubt about the suspension function of the Your® Sleeve, contact MotionTech Customer Service before allowing your patient to use the Your® Sleeve.

Trimming the product

You may trim the product to a shorter length. For a BK liner, you do not trim it shorter than 8cm (3.1") above the proximal end of the socket. For a sleeve, do not trim it shorter than 15cm (6") from mid-patella both towards proximal and distal extremities. The distal extremity of a Your® Sleeve from Scan should not be cut shorter, as this creates a risk of loss of suspension.

Use sharp scissors, a paper cutter, or a dedicated Hand-Held Cutting Wheel. To help prevent the seam from unraveling: Place a small drop of instant adhesive on the seam directly below the trimmed edge.

5. Miscellaneous

General Disclaimer

A Your® Liner/Sleeve can be ordered by a certified prosthetist only. A Your® Liner/Sleeve offers the possibility to the prosthetist to customise a wide variety of parameters, in order to provide optimal comfort and performances to his/her patient. It is the responsibility of the ordering prosthetist to select the suitable parameters depending on the patient s/he intends to provide a Your® Liner/Sleeve with.

The Your® Liner/Sleeve documentation gives advices on which parameters should be selected for particular situations. They

should be treated as guidelines, and not absolute rules, as each patient may differ. MotionTech's customer support teams are entirely at your disposal for further guidance on the Your® Liner/ Sleeve design process.

MotionTech retains all intellectual property rights reflected in its physical products, regardless of the transfer of the physical products to another party or parties.

Liability

Swiss Motion Technologies SA, as manufacturer, recommends using its products only under the specified conditions and for the intended purposes. A Your® Liner/Sleeve must be maintained according to the instructions for use. The manufacturer is not liable for any damage caused by combinations of component that were not specifically recommended by the manufacturer in the present document.

6. Ordering instructions

Important instructions for placing correct orders is provided in MotionTech's Frequently Asked Questions for professionals page (<https://www.motiontech.ch/faq-pros>).

MotionTech's ordering platform also include instructions and guidance on important steps to follow when placing orders (<https://www.motiontech.ch/yourlinerorderingplatform>).

Failure to properly follow these instructions may result in exposing your patients to risk of skin injuries such as blisters or other types of skin breakdown. Contact your local MotionTech representative for further guidance on proper processes for ordering a Your® Liner/Sleeve.

NOTE

The measurements taken for a Your® Sleeve must be taken with the joint flexed at same angle as the angle you ask us to produce the sleeve in. Failure to follow this predicament may generate a risk of suspension loss leading to fall, or other risks such a discomfort or skin injuries.

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