

Pediatrics / XXS sizes

Example 1



Example 2



Example 3



With Your® Liner

Limb



3D model



Your® Liner



Indication The residual limb exhibits an unusually small circumference.



Typical issues

No small-enough off-the-shelf liner available
Off-the-shelf liners come in limited sizes, and for some patients, it's simply impossible to find liners which are small-enough.
In such cases, upper-limb liners are sometimes used for lower limb patients, but such liners are not necessarily designed for weight bearing, and can therefore offer inadequate characteristics (e.g. too-thin, too-stiff material, suitable suspension options unavailable, etc.)



Key features



Anatomically tailored

Unlimited min size

Your® Liner does not per-se present minimum circumferences. Some design parameters (e.g. locking umbrella, thickness) have a lower boundary, but Your® Liner can still match even the smallest residual limbs, for which adequately sized off-the-shelf liners simply don't exist.

Min. umbrella diameter (pin-lock suspension): 90cm

Min. thickness: 2.5mm (Dynamic) ; 3mm (Active / Relax)

How can Your® Liner help ?



Added value: customisable aesthetics

The aesthetic aspect of prosthetic equipment can have an impact on how it is perceived and accepted by patients. The possibility of customising the aesthetics of the liner (and socket) offered by Your® Liner can be perceived as extra value to facilitate prosthesis acceptance.

How can Your® Liner help ?

